



**BANBROS ENGINEERING PVT. LTD.**

The Precision Measurement People

# METALLOGRAPHY SPECIMEN GRIDING & POLISHING MACHINE



- *VALUE*
- *VERSATILITY*
- *PERFORMANCE*

**BP-2G SERIES**

**DESCRIPTION:**

Model **BP-2G** is a cabinet type metallography specimen polishing machine. This machine features nice appearance design, big store cabinet and cooling system (which can clear up the heat produced during grinding so as to prevent the metallography structure from being damaged).

By equipping the motor with variable rotatingspeed, it can meet the different demands ofdifferent materials when polishing. It is the idealspecimen preparing instrument for using infactories, scientific research institutes and laboratories of colleges.

**TECHNICAL SPECIFICATIONS:**

Model No.	BP-2G	BP-2G(V)
No. of Polishing disc	Double	Double
Diameter of polishing disc	200mm	200mm
Rotating disk type	Fixed	Variable
Polishing disk rotating speed	900 rpm	50-900 rpm
Motor	2x80W 220V 50Hz.	
Net weight	43 Kgs	

\* Due to continues product development, product image & Specifications are subject to change.

\* Customization can be done as per user requirement.