



STEREO ZOOM MICROSCOPE

- *VALUE*
- *VERSATILITY*
- *PERFORMANCE*



BSZ-41/41T

**DESCRIPTION:**

Stereo Zoom Microscopes produces a positive image. They are adopted with horizontal knob for zoom control; they have long working distance, clear and smooth image quality, wide field, convenience operation etc. They can be used for culture and education, scientific research, farming and forestry, electronics and precision machine industries widely.

TECHNICAL SPECIFICATION:

Model No.	BSZ-41	BSZ-41T
Viewing Head	Binocular	Trinocular
Eyepieces	WF 10X (20mm)	
Eyepieces Tubes	45° and rotated 360°	
Objectives	1x / 4x	
Zoom Ratio	1:4	
Total Magnification	10x to 40x (Zoom)	
Field of View (mm)	23 to 5.5	
Working Distance(mm)	91	
Working Stage	Φ 95mm Frosted Glass	
Illuminator	Incident: 12V, 10W Transmitted: 12V, 10W	

* Due to continuous product development, Image & specification can be upgrade.