



BANBROS ENGINEERING PVT. LTD.

The Precision Measurement People

STEREO ZOOM MICROSCOPE



BSZ-608T

**FEATURES:**

- The leading role of the same level microscopes, wide range observation with zoom range 0.6X-5X,
- Zoom ratio 1:8.3 and visual field up to $\phi 76\text{mm}$.
- With high quality optical components, providing sharp and striking contrast image and no ghost image.
- Ease the tiredness with ergonomic design.
- H With Complete eye pieces, auxiliary objectives and other accessories, it is an ideal instrument for Industry assembling, inspection and teaching field.

TECHNICAL SPECIFICATIONS:

Model No.	BSZ-608T
Viewing Tube	30° Trinocular viewing head, Inclined at 30°
Diopter Adjustment Eyepiece	WF10X/ $\phi 23$ Extra wide field eyepiece WF10X/ $\phi 23$
Zoom Objective	0.6X-5X
Magnification	6X to 50X
Zoom Ratio	1:8.3
Working Distance	115mm
Interpupillary Distance	55-75mm
Stage Insert	Glass plate, Diameter 125mm
	Opal & black acrylic plate, Diameter 125mm
Focusing Range	105mm
Illumination	100V-240V/LED Incident Illumination
	100V-240V/LED Transmitted Illumination
Video Adapter	1x C Mount

* Due to continuous product development, Image & specification can be upgrade.