



METALLOGRAPHIC SPECIMEN GRINDING & POLISHING MACHINE

- *VALUE*
- *VERSATILITY*
- *PERFORMANCE*



BTS-8202



DESCRIPTION:

Model **MBP-2B** metallographic specimen grinding & polishing machine is equipped with double disc and features stepless speed changing grinding and polishing.

The rotating speed of grinding and polishing disc is adjustable from 50 to 1000 rpm. It can realize the whole process of rough grinding, fine grinding, rough polishing and finishing polishing for specimen preparing.

This machine features easy operation and economical applying. It is the ideal specimen preparing instrument for using in factories, scientific research institutes and laboratories of colleges.

TECHNICAL SPECIFICATIONS:

| Model No. | MBP-2BP |
|-------------------------------|---------------------|
| No. of Polishing disc | Double disc |
| Polishing head | With polishing head |
| Display of rotating speed | Digital display |
| Diameter of polishing disc | Ø 200mm |
| Diameter of grinding disc | Ø230mm |
| Polishing disc rotation | Variable |
| Polishing disc rotating speed | 50-1000 RPM |
| Motor | 550W 220V 50Hz |
| Dimensions | 690 X 715 X 310mm |
| Net weight | 45 Kg (approx.) |

* Image is for reference purpose only, due to continuous product development and upgradation the image and specification are subject to change.