



**BANBROS ENGINEERING PVT. LTD.**

The Precision Measurement People

# METALLOGRAPHIC SPECIMEN MOUNTING PRESS

- **VALUE**
- **VERSATILITY**
- **PERFORMANCE**



**ABQ2B-50**

**DESCRIPTION:**

**ABQ2B-50** This machine is Automatic Mounting Press for metallographic specimen. It is suitable for all thermosetting materials, because of the cooling system. It can automatically finish with the heating temperature, the holding time and acting force. The operator can change 4 specifications patterns according to the difficult specimens it can make 2 specimens at one time.

**TECHNICAL SPECIFICATIONS:**

Model No.	ABQ2B-50
Pattern Specimen	Ø25mm, Ø30mm, Ø40mm, Ø50mm
Heater Capacity	650W, 220V
Max. Power Consumption	1800W
Adjustable Temperature Range	0-190°C
Cooling Method	Water Cooling
Air Pressure	0-2MPa (corresponding preparing pressure intensity: 0-72MPa)
Power	220V, 50Hz, Single Phase
Heat Preserving Time	0-99mins and 99 Seconds
Dimension	615 × 510 × 500 cm
Weight	110Kg

\* Due to continues product development, image & Specifications are subject to change.